

Most Advanced Access System

- ▶ Advanced seal design
- ▶ Smooth instrument exchange
- ▶ Sequential Helix® technology

15mm Kii® Optical

Access System



15mm Kii® Optical Access System



Unparalleled Optical Clarity

Provides clear visualization of tissue layers without glare

Advanced Seal Design

Offers high flow insufflation and exceptional durability

Facilitates simple seal detachment for rapid desufflation



Optimized Cannula

Accommodates gastric bands, large staplers and 5mm instruments

Separator® Technology

Enables true non-bladed access with minimal fascial defect by separating muscle fibers along their natural planes



Funneled Entry

Ensures smooth instrument exchanges with minimal drag

REORDER #	DESCRIPTION	QUANTITY
C0R36	15mm x 100mm, Non-Threaded Kii Optical Access System	6/BOX
C0R37	15mm x 100mm, Threaded Kii Optical Access System	6/BOX
C0R38	15mm x 150mm, Non-Threaded Kii Optical Access System*	6/BOX
C0R39	15mm x 150mm, Threaded Kii Optical Access System*	6/BOX

* Please contact your Sales Representative for availability.

© 2009 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design, Separator and Kii are U.S. registered trademarks of Applied Medical Resources Corporation. SC01396A

Applied Medical
The New Generation Medical Company

Rancho Santa Margarita, CA 92688
Tel 949.713.8300
Fax 949.713.8400
www.appliedmedical.com

To place an order, or for more information, call 800.282.2212