

## Enabling Spaceguard™ Technology

- ▶ Minimizes depth of penetration necessary to initiate insufflation
- ▶ Establishes pneumoperitoneum under visualization with Fios obturator

# Kii® Fios®

First Entry



# Kii® Fios® First Entry

## Laparoscope Seal

*Maintains pneumoperitoneum with or without the use of a scope*


## Kii Seal Technology

*Allows smooth instrument exchanges, high flow insufflation and exceptional durability*



## Spaceguard™ Technology

*Provides rapid insufflation with minimal penetration into the peritoneal cavity*



*Insufflation apertures*



## Separator® Tip Technology

*Enables true non-bladed access with minimal fascial defect*

*Scope is not included*

REORDER NO.		DESCRIPTION	QUANTITY
NON-THREADED	THREADED		
CFI03	CFI13	5mm x 100mm Kii Fios First Entry	6/BOX
CFI04	CFI14	5mm x 100mm Kii Fios First Entry, Handled*	6/BOX
CFI33	CFI43	11mm x 100mm Kii Fios First Entry*	6/BOX
CFI34	CFI44	11mm x 100mm Kii Fios First Entry, Handled*	6/BOX
CFI73	CFI83	12mm x 100mm Kii Fios First Entry*	6/BOX
CFI74	CFI84	12mm x 100mm Kii Fios First Entry, Handled*	6/BOX

\* Please contact your Sales Representative for availability.

© 2008 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design, Fios, Kii, and Separator are registered trademarks of Applied Medical Resources Corporation.  
Spaceguard is a trademark of Applied Medical Resources Corporation. SC01383A

**Applied Medical**  
*The New Generation Medical Company*

Rancho Santa Margarita, CA 92688  
Tel 949.713.8300  
Fax 949.713.8400  
www.appliedmedical.com

To place an order, or for more information, call 800-282-2212