

## Precision and Strength

- ▶ Dual-edge blades
- ▶ Low profile mechanism
- ▶ Sustained sharpness over extended use

# Direct Drive®

Laparoscopic Scissors




# Direct Drive® Laparoscopic Scissors

## Precision and Strength




### Dual-Edge Blade

*Enables division of fibrotic tissue*  
*Allows transection and dissection of fine tissue*  
*Sustains sharpness over extensive use*



### Monopolar Electrosurgery Connector

*Utilizes free-swivel post sleeve to keep ESU cord off scissor handle*



### 360° Rotation Dial

*Allows convenient, single digit manipulation*

### Low Profile Mechanism

*Eliminates slots and protrusions, ensuring insulation down to the blades*

### Ergonomic Handle

*Increases comfort with soft finger inserts*

REORDER NO.	DESCRIPTION	QUANTITY
CB030	5mm x 35cm, DIRECT DRIVE DISPOSABLE LAPAROSCOPIC SCISSORS	10/BOX
CB040	5mm x 45cm, DIRECT DRIVE DISPOSABLE LAPAROSCOPIC SCISSORS	10/BOX



Rancho Santa Margarita, CA 92688  
Tel 949.713.8300  
Fax 949.713.8400  
www.appliedmedical.com